

Substanzname	Berichtsgrenze/ reporting limit (RL) [mg/kg]	Bentazon (Summenparameter)**	Chloroxuron	0,010	Dichlobenil	0,010	
		Bentazon-6-OH	0,010	Chlorpropham	0,010	Dichlofenthion	0,010
		Bentazon-8-OH	0,010	Chlorpropylat	0,010	Dichlofluanid	0,010
1,2,3,4-Tetrachlorbenzol	0,010	Benthiavalicarb		Chlorpyrifos	0,010	Dichlorvos	0,010
1,2,4,5-Tetrachlorbenzol	0,010	(Summenparameter)**		Chlorpyrifos-Methyl	0,010	Diclobutrazol	0,010
1,4-Dimethylnaphthalin	0,010	Benthiavalicarb-Isopropyl	0,010	Chlorsulfuron	0,010	Diclofop	0,010
1-Naphthylacetamid (1-NAD)	0,010	Benzovindiflupyr	0,010	Chlorthal-Dimethyl	0,010	Diclofop-Methyl	0,010
1-Naphthylelessigsäure	0,010	Benzoylprop-Ethyl	0,010	Chlorthalonil	0,010	Diclofop-Methyl	
1-Naphthylelessigsäure (Summenparameter)**		Benzyladenin	0,010	Chlorthiamid	0,010	(Summenparameter)**	
2,3,5-Trimethacarb/3,4,5- Trimethacarb	0,010	Bifenazat	0,010	Chlorthion	0,010	Dicloran	0,010
2,4-D (freie Säure)*	0,010	Bifenazat		Chlorthiophos	0,010	Dicofol (Summe aus p,p- und o,p-Isomeren)	0,010
2,4-D (Summenparameter)**		Bifenazat-diazen	0,010	Chlortoluron	0,010	Dicrotophos	0,010
2,4-DB (freie Säure)*	0,010	Bifenox	0,010	Chlozolinat	0,010	Dieldrin	0,010
2,4-DB (Summenparameter)**		Bifenthrin (Summe der Isomeren)	0,010	Cinerin I	0,010	Dieldrin (Summenparameter)**	
2,4-DP (freie Säure)*	0,010	Binapacryl	0,010	Cinerin II	0,010	Diethofencarb	0,010
2,4-DP (Summenparameter)**		Biphenyl	0,010	Cinidon-Ethyl (Cinidon-Ethyl und E-Isomer)	0,010	Diethyltoluamid (DEET)	0,010
3-Chloranilin	0,010	Bispyribac	0,010	Cinosulfuron	0,010	Difenoconazol	0,010
3-Decen-2-On	0,010	Bitertanol (Summe der Isomeren)	0,010	Clethodim	0,010	Difenoxyuron	0,010
479M04	0,010			Clethodim		Diflubenzuron	0,010
479M08	0,010	Bixafen	0,010	(Summenparameter)**		Diflufenican	0,010
479M16	0,010	Boscalid	0,010	Clethodim sulfon	0,010	Dimefox	0,010
4-CPA	0,010	Bromacil	0,010	Clethodim sulfoxid	0,010	Dimefuron	0,010
8-Hydroxychinolin	0,010	Bromfenvinphos	0,010	Climbazol	0,010	Dimethachlor	0,010
Abamectin (Summenparameter)**		Bromocyclen	0,010	Clodinafop		Dimethenamid (Summe aller Isomere)	0,010
Acephat	0,010	Bromophos-Ethyl	0,010	(Summenparameter)**			
Acequinocyl	0,010	Bromophos-Methyl	0,010	Clodinafop-Propargyl	0,010	Dimethoat	0,010
Acetamiprid	0,010	Bromoxynil	0,010	Clofentezin	0,010	Dimethomorph (Summe der Isomere)	0,010
Acetochlor	0,010	Brompropylat	0,010	Clomazon	0,010		
Acibenzolar-Säure (freie Säure)*	0,010	Bromuconazol	0,010	Clopyralid	0,010	Dimethylaminosulfotoluidid	0,010
Acibenzolar-S-methyl		(Diastereoisomere)		Cloquintocet-mexyl	0,010	Dimethylphenylformamide, 2,4-	
Acibenzolar-S-methyl (Summenparameter)**	0,010	BTS 44595	0,010	Clothianidin	0,010		
Acclonifen	0,010	BTS 44596	0,010	Coumaphos	0,010	Dimethylphenyl-N- methylformidin, N-2,4-	0,010
Acrinathrin	0,010	Bupirimat	0,010	Crimidin	0,010		
Alachlor	0,010	Buprofezin	0,010	Crotoxyphos	0,010	Dimetilan	0,010
Aldicarb	0,010	Butachlor	0,010	Crufomat	0,010	Dimoxystrobin	0,010
Aldicarb (Summenparameter)**		Butafenacil	0,010	Cyanazin	0,010	Diniconazol (Summe der Isomeren)	0,010
Aldicarb-sulfoxid	0,010	Butoxycarboxim	0,010	Cyanofenphos	0,010		
Aldoxycarb	0,010	Butralin	0,010	Cyanophos	0,010	Dinitramin	0,010
Aldrin	0,010	Buturon	0,010	Cyantraniliprol	0,010	Dinoseb	0,010
Allethrin	0,010	Cadusafos	0,010	Cyazofamid	0,010	Dinoseb (Summenparameter)**	
Ametoctradin	0,010	Captan	0,010	Cyclanilid	0,010	Dinoseb-Acetat	0,010
Ametryn	0,010	Captan (Summenparameter)**		Cycloat	0,010	Dinotefuran	0,010
Amidosulfuron	0,010	Carbaryl	0,010	Cyflufenamid (Summe der Isomeren)	0,010	Dioxacarb	0,010
Aminocarb	0,010	Carbendazim	0,010	Cyflumetofen	0,010	Dioxathion	0,010
Amisulbrom	0,010	Carbendazim		Cyfluthrin (Summe aller Isomeren)	0,010	Diphenamid	0,010
Amitraz	0,010	(Summenparameter)**				Diphenylamin	0,010
Amitraz (Summenparameter)**		Carbetamid (Summe der Isomeren)	0,010	Cyhalofop-butyl	0,010	Dipropetryn	0,010
Ancymidol	0,010	Carbofuran	0,001	Cyhalothrin, lambda-/gamma-	0,010	Disulfoton	0,010
Anilazin	0,010	Carbofuran		Cyhexatin (Summe Azocyclotin und Cyhexatin)	0,010	Disulfoton	
Atrazin	0,010	Carbofuran, 3-hydroxy	0,001	Cymiazol	0,010	(Summenparameter)**	
Atrazin, Desethyl-	0,010	Carbophenothion	0,010	Cymoxanil	0,010	Disulfoton-Sulfon	0,010
Avermectin B1a	0,006	Carbophenothion-Methyl	0,010	Cypermethrin (Summe aller Isomeren)	0,010	Disulfoton-Sulfoxid	0,010
Avermectin B1b	0,006	Carbosulfan	0,010	Cyproconazol	0,010	Ditalimfos	0,010
Azaconazol	0,010	Carboxin	0,010	Cyprodinil	0,010	Dithianon	0,010
Azadirachtin	0,010	Carboxin		Cyprofuram	0,010	Diuron	0,010
Azamethiphos	0,010	(Summenparameter)**		Cyromazin	0,010	DNOC	0,010
Azinphos-Ethyl	0,010	Carboxin-Sulfoxid	0,010	DDAC (Summenparameter)**		Dodemorph	0,010
Azinphos-Methyl	0,010	Carfentrazon-Ethyl	0,010	DDAC 10	0,010	Dodin	0,010
Aziprotryn	0,010	Cetrimoniumchlorid	0,010	DDAC 12	0,010	Edifenphos	0,010
Azoxystrobin	0,010	Chinomethionat	0,010	DDAC 8	0,010	Emamectinbenzoat B1a (als Emamectin)	0,002
BAC (Summenparameter)**		Chlorantraniliprol	0,010	DDD, p,p'-	0,010	Endosulfan	0,010
BAC 10	0,010	Chlorbensid	0,010	DDE, p,p'-	0,010	(Summenparameter)**	
BAC 12	0,010	Chlorbenzilat	0,010	DDT (Summenparameter)**		Endosulfan, alpha-	0,010
BAC 14	0,010	Chlorbromuron	0,010	DDT, o,p'-	0,010	Endosulfan, beta-	0,010
BAC 16	0,010	Chlorbufam	0,010	DDT, p,p'-	0,010	Endosulfansulfat	0,010
BAC 18	0,010	Chlordan (cis-/trans-)**		DEF	0,010	Endrin	0,010
BAC 8	0,010	Chlordan, cis-	0,010	Deltamethrin	0,010	EPN	0,010
Benalaxyl (Summe der Isomeren)	0,010	Chlordan, trans-	0,010	Demeton-S-Methyl	0,010	Epoxiconazol	0,010
Bendiocarb	0,010	Chlorfenapyr	0,010	Demeton-S-Methylsulfon	0,010	EPTC	0,010
Benfluralin	0,010	Chlorfenprop-Methyl	0,010	Denatoniumbenzoat	0,010	Etaconazol	0,010
Benfuracarb	0,010	Chlorfensulfon	0,010	Desethyl-Simazin	0,010	Ethiofencarb	0,010
Benodanil	0,010	Chlorfenvinphos	0,010	Desmedipham	0,010	Ethiofencarb-sulfon	0,010
Bensulfuron-Methyl	0,010	Chlorfluaazuron	0,010	Desmetryn	0,010	Ethion	0,010
Bentazon	0,010	Chloridazon	0,010	Diafenthiuron	0,010	Ethiprol	0,010
		Chlormephos	0,010	Diallat	0,010	Ethirimol	0,010
		Chloroneb	0,010	Diazinon	0,010	Ethofumesat	0,010
				Dicamba	0,010	Ethofumesat (Summenparameter)**	
						Ethofumesat-2-keto	0,010

Ethofumesat-Carbonsäure	0,010	Flufenacet thioglycolate	0,010	Isazofos	0,010	Mevinphos (Summe der E- und Z-Isomeren)	0,010
Ethoprophos	0,010	sulfoxid		Isobenzan	0,010	Milbemectin	
Ethoxyquin	0,010	Flufenoxuron	0,010	Isocarbamid	0,010	(Summenparameter)**	
Etofenprox	0,010	Flumethrin	0,010	Isocarbophos	0,010	Milbemectin A3	0,010
Etiozaxol	0,010	Flumetralin	0,010	Isodrin	0,010	Milbemectin A4	0,010
Etrifos	0,010	Flumioxazin	0,010	Isofenphos	0,010	Mirex	0,010
Famoxadon	0,010	Fluometuron	0,010	Isofenphos-Methyl	0,010	Molinat	0,010
Famphur	0,010	Fluopicolid	0,010	Isofetamid	0,010	Monocrotophos	0,010
Fenamidon	0,010	Fluopyram	0,010	Isomethiozin	0,010	Monolinuron	0,010
Fenamiphos	0,010	Fluotrimazol	0,010	Isopropcarb	0,010	Monuron	0,010
Fenamiphos		Fluoxastrobin	0,010	Isopropalin	0,010	Myclobutanil	0,010
(Summenparameter)**		Flupyradifuron	0,010	Isoprothiolan	0,010	Naled	0,010
Fenamiphos-sulfon	0,010	Fluquinconazol	0,010	Isoproturon	0,010	Napropamid	0,010
Fenamiphos-sulfoxid	0,010	Flurochloridon	0,010	Isopyrazam	0,010	Neburon	0,010
Fenarimol	0,010	Fluroxypyr	0,010	Isoxaben	0,010	Nicosulfuron	0,010
Fenazaquin	0,010	Fluroxypyr		Isoxadifen-Ethyl	0,010	Nitenpyram	0,010
Fenbuconazol	0,010	(Summenparameter)**		Isoxaflutol		Nitralin	0,010
Fenbutatinoxid	0,010	Flurprimidol	0,010	Isoxaflutol		Nitrapyrin	0,010
Fenchlorphos	0,010	Flusilazol	0,010	(Summenparameter)**		Nitrofen	0,010
Fenchlorphos		Fluthiacet-Methyl	0,010	Isoxaflutol diketonitril	0,010	Nitrothal-Isopropyl	0,010
(Summenparameter)**		Flutianil	0,010	Isoxathion	0,010	Norflurazon	0,010
Fenchlorphos-oxon	0,010	Flutolanil	0,010	Jasmolin I	0,010	Novaluron	0,010
Fenfluthrin	0,010	Flutriafol	0,010	Jasmolin II	0,010	Nuarimol	0,010
Fenhexamid	0,010	Fluvalinat (Summe der Isomere)	0,010	Karanjin	0,010	Ofurac	0,010
Fenitrothion	0,010	Fluxapyroxad	0,010	Kresoxim-Methyl	0,010	Omethoat	0,010
Fenobucarb	0,010	FM-6-1	0,010	Lenacil	0,010	Oxadiazon	0,010
Fenoprop (2,4,5-TP)	0,010	Folpet	0,010	Leptophos	0,010	Oxadixyl	0,010
Fenothiocarb	0,010	Folpet (Summenparameter)**		Linuron	0,010	Oxamyl	0,010
Fenoxaprop-Ethyl	0,010	Fonofos	0,010	Lufenuron	0,010	Oxathiapiprolin	0,010
Fenoxycarb	0,010	Forchlorfenuron	0,010	Malaoxon	0,010	Oxycarboxin	0,010
Fenpicoxamid	0,010	Formetanat	0,010	Malathion	0,010	Oxychlordan	0,010
Fenpropathrin	0,010	Formetanat (als Formetanat-HCl)**		Malathion		Oxydemeton-Methyl	0,010
Fenpropidin	0,010	Formetanat		(Summenparameter)**		Oxydemeton-Methyl (Summenparameter)**	
Fenpropimorph (Summe der Isomeren)	0,010	Formothion	0,010	Mandestrobilin	0,010	Oxyfluorfen	0,010
Fenpyrazamin	0,010	Fosthiazat	0,010	Mandipropamid (Summe der Isomere)	0,010	Oxymatrin	0,010
Fenpyroximat	0,010	Fuberidazol	0,010	Matrin	0,010	Paclobutrazol (Summe der Isomeren)	0,010
Fenson	0,010	Furalaxyl	0,010	MCPA (freie Säure)*	0,010	Paraoxon	0,010
Fensulfothion	0,010	Furathiocarb	0,010	MCPA (Summenparameter)**		Parathion	0,010
Fensulfothion		Furmecyclox	0,010	MCPB (freie Säure)*	0,010	Parathion-Methyl	0,010
(Summenparameter)**		Genite	0,010	Mecarbam	0,010	(Summenparameter)**	
Fensulfothion-sulfon	0,010	Halfenprox	0,010	Mecoprop (Summe der Isomeren)	0,010	Pebulat	0,010
Fenthion	0,010	Haloxypfop (freie Säure)*	0,010	Mefenpyr-diethyl	0,010	Penconazol (Summe der Isomeren)	0,010
Fenthion		Haloxypfop		Mefentrifluconazol	0,010	Pencycuron	0,010
(Summenparameter)**		(Summenparameter)**		Mepanipyrim	0,010	Pencycuron	
Fenthion-Oxon	0,010	Haloxypfop-ethoxyethylster*	0,010	Mephosfolan	0,010	(Summenparameter)**	
Fenthion-Oxonsulfon	0,010	Haloxypfop-methylester*	0,010	Merphos	0,010	Pencycuron-PB-amin	0,010
Fenthion-Oxonsulfoxid	0,010	HCH, alpha-	0,010	Metaflumizon (E- und Z-Isomere)	0,010	Pendimethalin	0,010
Fenthion-sulfon	0,010	HCH, beta-	0,010	Metaxyl (Summe der Isomeren)	0,010	Penflufen	0,010
Fentin	0,010	HCH, delta-	0,010	Metaldehyd	0,010	Pentachloranilin	0,010
Fenuron	0,010	HCH, epsilon-	0,010	Metamitron	0,010	Pentachloranisol	0,010
Fenvalerat/Esfenvalerat (Summe)	0,010	HCH, gamma-	0,010	Metazachlor		Pentachlorbenzol	0,010
Fipronil	0,005	Heptachlor	0,010	(Summenparameter)**		Pentachlorphenol	0,010
Fipronil (Summenparameter)**		Heptachlor		Metconazol (Summe der Isomeren)	0,010	Pentanochlor	0,010
Fipronil-Sulfid	0,005	(Summenparameter)**		Methabenzthiazuron	0,010	Penthiopyrad	0,010
Fipronil-sulfon	0,005	Heptachlorepoxyd, cis-	0,010	Methacrifos	0,010	Permethrin (Summe der Isomeren)	0,010
Flamprop-M-Isopropyl	0,010	Heptachlorepoxyd, trans-	0,010	Methamidophos	0,010	Perthan	0,010
Flamprop-M-Methyl	0,010	Heptenophos	0,010	Methidathion	0,010	Pethoxamid	0,010
Flazasulfuron	0,010	Hexachlorbenzol	0,010	Methiocarb	0,010	Phenmedipham	0,010
Flonicamid	0,010	Hexaconazol	0,010	(Summenparameter)**		Phenothrin	0,010
Flonicamid		Hexaflumuron	0,010	Methiocarb-sulfon	0,010	Phenthoat	0,010
(Summenparameter)**		Hexazinon	0,010	Methiocarb-sulfoxid	0,010	Phenylphenol, ortho-	0,010
Florpyrauxifen-benzyl	0,010	Hexythiazox	0,010	Methomyl	0,010	Phorat (Summenparameter)**	
Fluacrypyrim	0,010	Hydramethylnon	0,010	Methoprotryn	0,010	Phorat-Oxon	0,010
Fluazifop (freie Säure)*	0,010	Icaridin	0,010	Methoxychlor	0,010	Phorat-Oxon-Sulfon	0,010
Fluazifop		Imazalil	0,010	Methoxyfenozid	0,010	Phorat-sulfon	0,010
(Summenparameter)**		Imazamox	0,010	Metobromuron	0,010	Phosalon	0,010
Fluazifop-Butyl*	0,010	Imazapyr	0,010	Metolachlor (Summe der Isomeren)	0,010	Phosfolan	0,010
Fluazinam	0,010	Imazaquin	0,010	Metolcarb	0,010	Phosmet	0,010
Fluazuron	0,010	Imazethapyr	0,010	Metominostrobin	0,010	Phosmet (Summenparameter)**	
Flubendiamid	0,010	Imibenconazol	0,010	Metoxuron	0,010	Phosmet Oxon	0,010
Fluchloralin	0,010	Imidacloprid	0,010	Metrafenon	0,010	Phosphamidon	0,010
Flucythrinate	0,010	Indoxacarb (Summe der R- und S-Isomeren)	0,010	Metribuzin	0,010	Phoxim	0,010
Fludioxonil	0,010	Iodofenphos	0,010	Metsulfuron-Methyl	0,010	Phthalimid	0,010
Flusulfuron	0,010	Iodosulfuron-Methyl	0,010			Picloram	0,010
Flufenacet	0,010	Ioxynil	0,010				
Flufenacet		Ioxynil (Summenparameter)**					
(Summenparameter)**		Iprobenfos	0,010				
Flufenacet oxalat	0,010	Iprodion	0,010				
Flufenacet Sulfonsäure	0,010	Iprovalicarb	0,010				

Picolinafen	0,010	Saflufenacil, M800H35	0,010	Tricyclazol	0,010	Chlormequat, Mepiquat
Picoxystrobin	0,010	Sebutylazin	0,010	Tridemorph	0,010	LC-MS/MS
Piperonylbutoxid	0,010	Sethoxydim	0,010	Trietazin	0,010	SOP-Nr. 495:2016-10; EU-SRM QuPpe 2019-05
Piperophos	0,010	Silaneophan	0,010	Trifloxystrobin	0,010	BG/RL je 0,01 mg/kg
Pirimicarb	0,010	Silthiofam	0,010	Triflurizol	0,010	Chlormequat, BG/RL 0,01 mg/kg
Pirimiphos-Ethyl	0,010	Simazin	0,010	(Summenparameter)**		Mepiquat, BG/RL 0,01 mg/kg
Pirimiphos-Methyl	0,010	Sintofen	0,010	Triflurizol		Chlormequat (Summenparameter)**
Plifenat	0,010	Spinetoram	0,010	Triflurizol		Mepiquat (Summenparameter)**
Prochloraz	0,010	Spinosad		Trifluralin	0,010	Dithiocarbamate, berechnet als CS ₂ /calculated as CS ₂
Prochloraz		(Summenparameter)**		Triflursulfuron-Methyl	0,010	gem. ASU §64 LFGB L 00.00-49/2 modif.: GC-MSD
(Summenparameter)**		Spinosyn A	0,010	Triforin	0,010	BG/RL 0,01 mg/kg,
Procymidon	0,010	Spinosyn D	0,010	Trinexapac	0,010	Ethephon BNN
Profenofos	0,010	Spirodiclofen	0,010	Triticonazol	0,010	LC-MS/MS
Profluralin	0,010	Spiromesifen	0,010	Tritosulfuron	0,010	SOP-Nr. 495:2016-10; EU-SRM QuPpe 2019-05
Profoxydim	0,010	Spirotetramat	0,010	Uniconazol	0,010	BG/RL 0,01 mg/kg
Promecarb	0,010	Spirotetramat		Valifenalat	0,010	Ethylenoxid (Summenparameter)**
Prometon	0,010	(Summenparameter)**		Vamidothion	0,010	GC-MSMS, SOP-0653:02-2022
Prometryn	0,010	Spirotetramat-enol	0,010	Vinclozolin	0,010	BG/RL 0,01mg/kg
Propachlor	0,010	Spiroxamin (Summe der Isomeren)	0,010	Zoxamid	0,010	Ethylenoxid, 2-Chlorethanol, Ethylenoxid (Summe)
Propachlor		Sulfentrazon	0,010			GC-MSMS, SOP-0653:02-2022
(Summenparameter)**		Sulfosulfuron	0,010			BG/RL 0,01mg/kg
Propamocarb	0,010	Sulfotep	0,010			Ethylenoxid (Summenparameter)**
Propanil	0,010	Sulfoxaflor	0,010			Ethylenoxid, BG/RL 0,01 mg/kg
Propaquizafop*	0,010	Sulprofos	0,010			2 Chlorethanol, BG/RL 0,01 mg/kg
Propargit	0,010	Swep	0,010			Fosetyl (Summenparameter)**
Propazin	0,010	Tebuconazol	0,010			LC-MS/MS
Propetamphos	0,010	Tebufenozid	0,010			SOP-Nr. 495:2016-10; EU-SRM QuPpe 2019-05
Propham	0,010	Tebufenpyrad	0,010			Fosetyl, BG/RL 0,01 mg/kg
Propiconazol	0,010	Tebupirimphos	0,010			Phosphonsäure, BG/RL 0,01 mg/kg
Propoxur	0,005	Tebutam	0,010			Gesamtbromid
Propyzamid	0,010	Tecnazen	0,010			GC-MSD
Proquinazid	0,010	Teflubenzuron	0,010			ASU L 00.00-36/1:2004-07,
Prosulfocarb	0,010	Tefluthrin (Summe der Isomere)	0,010			BG/RL 1 mg/kg
Prosulfuron	0,010	Tembotriol	0,010			Gesamtbromid, Chlorid, Br/Cl-Verhältnis
Prothioconazol	0,010	Temephos	0,010			GC-MSD
Prothioconazol-Desthio	0,010	Tepraloxymid	0,010			ASU L 00.00-36/1:2004-07,
Prothioconazol-Desthio		Terbacil	0,010			Gesamtbromid, BG/RL 1 mg/kg
(Summenparameter)**		Terbufos	0,010			Chlorid 5 mg/kg
Prothiophos	0,010	Terbutryn	0,010			Bromid-Chlorid-Verhältnis
Pymetrozin	0,010	Tetrachlorvinphos	0,010			Glyphosat, AMPA, Glufosinat
Pyraclostrobin	0,010	Tetraconazol	0,010			LC-MS/MS
Pyraflufen	0,010	Tetradifon	0,010			SOP-0657:2022-12
Pyraflufen-Ethyl	0,010	Tetrahydrophthalimid	0,010			BG/RL je 0,01 mg/kg
Pyraflufen-Ethyl		Tetramethrin	0,010			Nikotin
(Summenparameter)**		Tetrasul	0,010			LC-MS/MS
Pyrazophos	0,010	TFNA	0,010			SOP-196: 2018-07
Pyrethrin		TFNG	0,010			BG/RL 0,01mg/kg
(Summenparameter)**		Thiabendazol	0,010			Nitrat
Pyrethrin I	0,010	Thiacloprid	0,010			HPLC/IC
Pyrethrin II	0,010	Thiamethoxam	0,010			Gemäß SOP-0570, Version 2.0; June 2020
Pyridaben	0,010	Thifensulfuron-Methyl	0,010			BG/RL 1 mg/kg
Pyridafol (CL 9673)	0,010	Thiodicarb	0,010			Paraquat, Diquat
Pyridalyl	0,010	Thiofanox	0,010			mittels LC-MS/MS
Pyridaphenthion	0,010	Thiofanox-sulfon	0,010			BG/RL je 0,01 mg/kg
Pyridat	0,010	Thiofanox-sulfoxid	0,010			Perchlorat, Chlorat
Pyrifenoxy	0,010	Thiometon	0,010			LC-MS/MS
Pyrimethanil	0,010	Thionazin	0,010			SOP-Nr. 495:2016-10; EU-SRM QuPpe 2019-05
Pyriofenon	0,010	Thiophanat-Methyl	0,010			BG/RL je 0,01 mg/kg
Pyriproxifen	0,010	Tiocardbazil	0,010			Phosphan
Quinalphos	0,010	Tolclofos-Methyl	0,010			GC-MSD
Quinclorac	0,010	Tolfenpyrad	0,010			SOP-0559:2019-05
Quinmerac	0,010	Tolylfluand	0,010			BG/RL 0,01 mg/kg
Quinoclam	0,010	(Summenparameter)**				Pyridat, Pyridafol, Pyridat (Summenparameter)**
Quinoxifen	0,010	Transfluthrin	0,010			nach Hydrolyse
Quintozen	0,010	Triadimefon	0,010			LC-MS/MS
(Summenparameter)**		Triadimenol (Summe der Isomeren)	0,010			DIN EN 15662:2018-07 Modulares QuEChERS-
Quizalofop (Summe der Isomeren) (freie Säure)*	0,010	Triallat	0,010			Verfahren, alkalische Hydrolyse
Quizalofop		Triamphos	0,010			BG/RL 0,01mg/kg
(Summenparameter)**		Triasulfuron	0,010			Sulfat
Quizalofop-Ethyl*	0,010	Triazamat	0,010			Gem. § 64 LFGB ASU L 00.00 46/1, November 1999
Rabenzazol	0,010	Triazophos	0,010			BG/RL 10 mg/kg
Resmethrin (Summe der Isomere)	0,010	Tribenuron-Methyl	0,010			Trimethylsulfonium (Trimesium)
Rimsulfuron	0,010	Trichlorfon	0,010			LC-MS/MS
Rotenon	0,010	Trichloronat	0,010			SOP-Nr. 495:2016-10; EU-SRM QuPpe 2019-05
Saflufenacil	0,010	Triclopyr	0,010			BG/RL 0,01mg/kg
Saflufenacil						
(Summenparameter)**						
Saflufenacil, M800H11	0,010					