

Substanzname	Berichtsgrenze/ reporting limit (RL) [mg/kg]						
		Benthiavalicarb-Isopropyl	0,010	Chlorthion	0,010	Diethofencarb	0,010
		Benzovindiflupyr	0,010	Chlorthiophos	0,010	Diethyltoluamid (DEET)	0,010
		Benzoylprop-Ethyl	0,010	Chlortoluron	0,010	Difenoconazol	0,010
		Benzyladenin	0,010	Chlozolinat	0,010	Difenoxuron	0,010
1,2,3,4-Tetrachlorbenzol	0,010	BH 518-2	0,010	Cinerin I	0,010	Diflubenzuron	0,010
1,2,4,5-Tetrachlorbenzol	0,010	BH 518-4	0,010	Cinerin II	0,010	Diflufenican	0,010
1,4-Dimethylnaphthalin	0,010	Bifenazat	0,010	Cinidon-Ethyl (Cinidon-Ethyl und E-Isomer)	0,010	Dimefox	0,010
1-Naphthylacetamid (1-NAD)	0,010	Bifenazat (Summenparameter)**		Cinosulfuron	0,010	Dimefuron	0,010
1-Naphthyllessigsäure	0,010	Bifenazat-diazen	0,010	Clethodim	0,010	Dimethachlor	0,010
1-Naphthyllessigsäure (Summenparameter)****		Bifenox	0,010	Clethodim (Summenparameter)**		Dimethenamid (Summe aller Isomere)	0,010
2,3,5-Trimethacarb/3,4,5-Trimethacarb	0,010	Bifenthrin (Summe der Isomere)	0,010	Clethodim sulfon	0,010	Dimethoat	0,010
2,4-D (freie Säure)*	0,010	Binapacryl	0,010	Clethodim sulfoxid	0,010	Dimethomorph (Summe der Isomere)	0,010
2,4-D (Summenparameter)**		Biphenyl	0,010	Climbazol	0,010	Dimethylaminosulfotolidid	0,010
2,4-DB (freie Säure)*	0,010	Bispyribac (freie Säure)*	0,010	Clodinafop (Summenparameter)**		Dimethylphenylformamide, 2,4-	0,010
2,4-DB (Summenparameter)**		Bispyribac (Summenparameter)**		Clodinafop-Propargyl	0,010	Dimethylphenyl-N-methylformidin, N-2,4-	0,010
2,4-DP (freie Säure)*	0,010	Bitertanol (Summe der Isomere)	0,010	Clofentezin	0,010	Dimetilan	0,010
2,4-DP (Summenparameter)**		Bixafen	0,010	Clomazon	0,010	Dimoxystrobin	0,010
3-Chloranilin	0,010	Boscalid	0,010	Clopyralid	0,010	Diniconazol (Summe der Isomere)	0,010
3-Decen-2-On	0,010	Bromacil	0,010	Cloquintocet-mexyl	0,010	Dinitramin	0,010
479M04	0,010	Bromfenvinphos	0,010	Clothianidin	0,010	Dinoseb	0,010
479M08	0,010	Bromocyclen	0,010	Coumaphos	0,010	Dinoseb (Summenparameter)**	
479M16	0,010	Bromophos-Ethyl	0,010	Crimidin	0,010	Dinoseb-Acetat	0,010
4-Bromphenylharnstoff	0,010	Bromoxynil	0,010	Crotoxyphos	0,010	Dinotefuran	0,010
4-CPA	0,010	Brompropylat	0,010	Cruformat	0,010	Dioxacarb	0,010
8-Hydroxychinolin	0,010	Bromuconazol (Diastereoisomere)	0,010	Cyanazin	0,010	Dioxathion	0,010
Abamectin (Summenparameter)**		BTS 44595	0,010	Cyanofenphos	0,010	Diphenamid	0,010
Acephat	0,010	BTS 44596	0,010	Cyanophos	0,010	Diphenylamin	0,010
Acequinocyl	0,010	Bupirimat	0,010	Cyantraniliprol	0,010	Dipropetryn	0,010
Acetamidrid	0,010	Buprofezin	0,010	Cyazofamid	0,010	Disulfoton	0,010
Acetochlor	0,010	Butachlor	0,010	Cyclanilid	0,010	Disulfoton (Summenparameter)**	
Acibenzolar-Säure (freie Säure)*	0,010	Butafenacil	0,010	Cycloat	0,010	Disulfoton-Sulfon	0,010
Acibenzolar-S-methyl	0,010	Butoxycarboxim	0,010	Cyflufenamid (Summe der Isomere)	0,010	Disulfoton-Sulfoxid	0,010
Acibenzolar-S-methyl (Summenparameter)**		Butralin	0,010	Cyflumetofen (Summe der Isomere)	0,010	Ditalimfos	0,010
Aclonifen	0,010	Buturon	0,010	Cyfluthrin (Summe aller Isomere)	0,010	Dithianon	0,010
Acrinathrin	0,010	Cadusafos	0,010	Cyhalofop-butyl	0,010	Diuron	0,010
Alachlor	0,010	Captan	0,010	Cyhalothrin, lambda-/gamma-	0,010	DNOC	0,010
Aldicarb	0,010	Captan (Summenparameter)**		Cyhexatin (Summe Azocyclotin und Cyhexatin)	0,010	Dodemorph	0,010
Aldicarb (Summenparameter)**		Carbaryl	0,010	Cymiazol	0,010	Dodin	0,010
Aldicarb-sulfoxid	0,010	Carbendazim	0,010	Cymoxanil	0,010	Edifenphos	0,010
Aldoxycarb	0,010	Carbendazim (Summenparameter)**		Cypermethrin (Summe aller Isomere)	0,010	Emamectin B1a (freie Base)	0,002
Aldrin	0,010	Carbetamid (Summe der Isomere)	0,010	Cyproconazol	0,010	Endosulfan (Summenparameter)**	
Allethrin	0,010	Carbofuran	0,001	Cyprodinil	0,010	Endosulfan, alpha-	0,010
Ametoctradin	0,010	Carbofuran (Summenparameter)**		Cyprofuram	0,010	Endosulfan, beta-	0,010
Ametryn	0,010	Carbofuran, 3-hydroxy	0,001	Cyromazin	0,010	Endosulfansulfat	0,010
Amidosulfuron	0,010	Carbophenothion	0,010	DDAC (Summenparameter)**		Endrin	0,010
Aminocarb	0,010	Carbophenothion-Methyl	0,010	DDAC 10	0,010	EPN	0,010
Amisulbrom	0,010	Carbosulfan	0,010	DDAC 12	0,010	Epoxiconazol	0,010
Amitraz	0,010	Carboxin	0,010	DDAC 8	0,010	EPTC	0,010
Amitraz (Summenparameter)**		Carboxin (Summenparameter)**		DDD, p,p'-	0,010	Etaconazol	0,010
Ancymidol	0,010	Carboxin-Sulfoxid	0,010	DDE, p,p'-	0,010	Ethiofencarb	0,010
Anilazin	0,010	Carfentrazon-Ethyl	0,010	DDT (Summenparameter)**		Ethiofencarb-sulfon	0,010
Anthrachinon	0,010	Cetrimoniumchlorid	0,010	DDT, o,p'-	0,010	Ethion	0,010
Atrazin	0,010	Chinomethionat	0,010	DDT, p,p'-	0,010	Ethiprol	0,010
Atrazin, Desethyl-	0,010	Chlorantraniliprol	0,010	DEF	0,010	Ethirimol	0,010
Avermectin B1a	0,006	Chlorbensid	0,010	Deltamethrin	0,010	Ethofumesat	0,010
Avermectin B1b	0,006	Chlorbenzilat	0,010	Demeton-S-Methyl	0,010	Ethofumesat (Summenparameter)**	
Azaconazol	0,010	Chlorbromuron	0,010	Demeton-S-Methylsulfon	0,010	Ethofumesat-2-keto	0,010
Azadirachtin	0,010	Chlorbufam	0,010	Denatoniumbenzoat	0,010	Ethofumesat-Carbonsäure	0,010
Azamethiphos	0,010	Chlordan (cis-/trans-)**		Desethyl-Simazin	0,010	Ethoprophos	0,010
Azinphos-Ethyl	0,010	Chlordan, cis-	0,010	Desmedipham	0,010	Ethoxyquin	0,010
Azinphos-Methyl	0,010	Chlordan, trans-	0,010	Desmetryn	0,010	Etofenprox	0,010
Aziprotryn	0,010	Chlordecon	0,010	Diafenthiuron	0,010	Etoazol	0,010
Azoxystrobin	0,010	Chlorfenapyr	0,010	Diallat	0,010	Etrimfos	0,010
BAC (Summenparameter)**		Chlorfenprop-Methyl	0,010	Diazinon	0,010	Famoxadon	0,010
BAC 10	0,010	Chlorfenson	0,010	Dicamba	0,010	Famphur	0,010
BAC 12	0,010	Chlorfenvinphos	0,010	Dichlobenil	0,010	Fenamidon	0,010
BAC 14	0,010	Chlorfluazuron	0,010	Dichlofenthiion	0,010	Fenamiphos	0,010
BAC 16	0,010	Chloridazon	0,010	Dichlofluanid	0,010	Fenamiphos (Summenparameter)**	
BAC 18	0,010	Chloridazon (Summenparameter)**		Dichlorvos	0,010	Fenamiphos-sulfon	0,010
BAC 8	0,010	Chloridazon-Desphenyl	0,010	Diclobutrazol	0,010	Fenamiphos-sulfoxid	0,010
Benalaxyl (Summe der Isomere)	0,010	Chlormephos	0,010	Diclofop	0,010	Fenarimol	0,010
Bendiocarb	0,010	Chloroneb	0,010	Diclofop-Methyl	0,010	Fenazaquin	0,010
Benfluralin	0,010	Chloroxuron	0,010	Diclofop-Methyl (Summenparameter)**		Fenbuconazol	0,010
Benfuracarb	0,010	Chlorpropham	0,010	Dicloran	0,010	Fenbutatinoxid	0,010
Benodanil	0,010	Chlorpropylat	0,010	Dicofol (Summe aus p,p- und o,p-Isomere)	0,010	Fenchlorphos	0,010
Bensulfuron-Methyl	0,010	Chlorpyrifos	0,010	Dicrotophos	0,010	Fenchlorphos-oxon	0,010
Bentazon	0,010	Chlorpyrifos-Methyl	0,010	Dieldrin	0,010	Fenfluthrin	0,010
Bentazon (Summenparameter)**		Chlorsulfuron	0,010	Dieldrin (Summenparameter)**		Fenhexamid	0,010
Bentazon-6-OH*	0,010	Chlorthal-Dimethyl	0,010			Fenitrothion	0,010
Bentazon-8-OH*	0,010	Chlorthalonil	0,010				
Benthiavalicarb (Summenparameter)**		Chlorthiamid	0,010				

Fenobucarb	0,010	Furmecyclox	0,010	Mepanipyrim	0,010	Pentachloranisol	0,010
Fenoprop (2,4,5-TP)	0,010	Genite	0,010	Mephosfolan	0,010	Pentachlorbenzol	0,010
Fenothiocarb	0,010	Halfenprox	0,010	Mepronil	0,010	Pentachlorphenol	0,010
Fenoxaprop-Ethyl	0,010	Halosulfuron-Methyl	0,010	Merphos	0,010	Pentachlor	0,010
Fenoxycarb	0,010	Haloxypfop (freie Säure)*	0,010	Metaflumizon (E- und Z-Isomere)	0,010	Penthiopyrad	0,010
Fenpicoxamid	0,010	Haloxypfop (Summenparameter)**		Metalaxyl (Summe der Isomere)	0,010	Permethrin (Summe der Isomere)	0,010
Fenpropathrin	0,010	Haloxypfop-ethoxyethylester*	0,010	Metaldehyd	0,010	Perthan	0,010
Fenpropidin	0,010	Haloxypfop-methylester*	0,010	Metamitron	0,010	Pethoxamid	0,010
Fenpropimorph (Summe der Isomere)	0,010	HCH, alpha-	0,010	Metazachlor		Phenmedipham	0,010
Fenpyrazamin	0,010	HCH, beta-	0,010	(Summenparameter)**		Phenothrin	0,010
Fenpyroximat	0,010	HCH, delta-	0,010	Metconazol (Summe der Isomere)	0,010	Phenthoat	0,010
Fenson	0,010	HCH, epsilon-	0,010	Methabenzthiazuron	0,010	Phenylphenol, ortho-	0,010
Fensulfthion	0,010	HCH, gamma-	0,010	Methacrifos	0,010	Phenylphenol, ortho-	
Fensulfthion-sulfon	0,010	Heptachlor	0,010	Methamidophos	0,010	(Summenparameter)**	
Fenthion	0,010	Heptachlor (Summenparameter)**		Methidathion	0,010	Phorat	0,010
Fenthion (Summenparameter)**		Heptachlorepoxyd, cis-	0,010	Methiocarb	0,010	Phorat (Summenparameter)**	
Fenthion-Oxon	0,010	Heptachlorepoxyd, trans-	0,010	Methiocarb (Summenparameter)**		Phorat-Oxon	0,010
Fenthion-Oxonsulfon	0,010	Heptenophos	0,010	Methiocarb-sulfon	0,010	Phorat-Oxon-Sulfon	0,010
Fenthion-Oxonsulfoxid	0,010	Hexachlorbenzol	0,010	Methiocarb-sulfoxid	0,010	Phorat-sulfon	0,010
Fenthion-sulfon	0,010	Hexaconazol	0,010	Methomyl	0,010	Phosalon	0,010
Fentin	0,010	Hexaflumuron	0,010	Methoprotryn	0,010	Phosfolan	0,010
Fenuron	0,010	Hexazinon	0,010	Methoxychlor	0,010	Phosmet	0,005
Fenvalerat/Esfenvalerat (Summe)	0,010	Hexythiazox	0,010	Methoxyfenozid	0,010	Phosphamidon	0,010
Fipronil	0,005	Hydramethylnon	0,010	Metobromuron	0,010	Phoxim	0,010
Fipronil (Summenparameter)**		Hymexazol	0,010	Metobromuron		Phthalimid	0,010
Fipronil-Sulfid	0,005	Icaridin	0,010	(Summenparameter)**		Picloram	0,010
Fipronil-sulfon	0,005	Imazalil	0,010	Metolachlor (Summe der Isomere)	0,010	Picolinafen	0,010
Flazasulfuron	0,010	Imazamox	0,010	Metolcarb	0,010	Picoxystrobin	0,010
Flometoquin	0,010	Imazapyr	0,010	Metominostrobin	0,010	Piperonylbutoxid	0,010
Flonicamid	0,010	Imazaquin	0,010	Metoxuron	0,010	Piperophos	0,010
Flonicamid (Summenparameter)**		Imazethapyr	0,010	Metrafenon	0,010	Pirimicarb	0,010
Florpyrauxifen-benzyl	0,010	Imibenconazol	0,010	Metribuzin	0,010	Pirimiphos-Ethyl	0,010
Fluacrypyrim	0,010	Imidacloprid	0,010	Metsulfuron-Methyl	0,010	Pirimiphos-Methyl	0,010
Fluazifop (freie Säure)*	0,010	Indaziflam	0,010	Mevinphos (Summe der E- und Z-Isomere)	0,010	Plifenat	0,010
Fluazifop (Summenparameter)**		Indoxacarb (Summe der Isomere)	0,010	Milbemectin		Prallethrin	0,010
Fluazifop-Butyl*	0,010	Iodofenphos	0,010	(Summenparameter)**		Pretilachlor	0,010
Fluazinam	0,010	Iodosulfuron-Methyl	0,010	Milbemectin A3	0,010	Prochloraz	0,010
Fluazuron	0,010	Ioxynil	0,010	Milbemectin A4	0,010	Prochloraz (Summenparameter)**	
Flubendiamid	0,010	Iprobenfos	0,010	Mirex	0,010	Procymidon	0,010
Fluchloralin	0,010	Iprodion	0,010	Molinal	0,010	Profenofos	0,010
Flucythrinat	0,010	Iprovalicarb	0,010	Monocrotophos	0,010	Profluralin	0,010
Fludioxonil	0,010	Isazofos	0,010	Monolinuron	0,010	Profoxydim	0,010
Fluensulfon	0,010	Isobenzan	0,010	Monuron	0,010	Promecarb	0,010
Flufenacet	0,010	Isocarbamid	0,010	Myclobutanil	0,010	Prometon	0,010
Flufenacet (Summenparameter)**		Isocarbophos	0,010	Naled	0,010	Prometryn	0,010
Flufenacet oxalat	0,010	Isodrin	0,010	Napropamid	0,010	Propachlor	0,010
Flufenacet Sulfonsäure	0,010	Isofenphos	0,010	Neburon	0,010	Propachlor (Summenparameter)**	
Flufenacet thioglycolate sulfoxid	0,010	Isofenphos-Methyl	0,010	Nicosulfuron	0,010	Propamocarb	0,010
Flufenoxuron	0,010	Isofetamid	0,010	Nitenpyram	0,010	Propanil	0,010
Flumethrin	0,010	Isomethiozin	0,010	Nitralin	0,010	Propaquizafop	0,010
Flumetralin	0,010	Isoprocab	0,010	Nitralin	0,010	Propargit	0,010
Flumioxazin	0,010	Isopropalin	0,010	Nitrofen	0,010	Propazin	0,010
Fluometuron	0,005	Isoprothiolan	0,010	Nitrothal-Isopropyl	0,010	Propetamphos	0,010
Fluopicolid	0,010	Isoproturon	0,010	Norflurazon	0,010	Propham	0,010
Fluopyram	0,010	Isopyrazam	0,010	Novaluron (Summe der Isomere)	0,010	Propiconazol	0,010
Fluotrimazol	0,010	Isoxaben	0,010	Nuarimol	0,010	Propoxur	0,005
Fluoxastrobin	0,010	Isoxadifen-Ethyl	0,010	Ofurac	0,010	Propyzamid	0,010
Flupyradifuron	0,010	Isoxaflutol	0,010	Omethoat	0,010	Proquinazid	0,010
Fluquinconazol	0,010	Isoxaflutol (Summenparameter)**		Oxadiazon	0,010	Prosulfocarb	0,010
Flurochloridon	0,010	Isoxaflutol diketonitril	0,010	Oxadixyl	0,010	Prosulfuron	0,010
Fluroxypyr (freie Säure)*	0,010	Isoxathion	0,010	Oxamyl	0,001	Prothioconazol	0,010
Fluroxypyr (Summenparameter)**		Jasmodin I	0,010	Oxathiapirolin	0,010	Prothioconazol-Desthio	0,010
Flurprimidol	0,010	Jasmodin II	0,010	Oxycarboxin	0,010	Prothioconazol-Desthio	
Flusilazol	0,010	Karanjin	0,010	Oxychlordan	0,010	(Summenparameter)**	
Fluthiacet-Methyl	0,010	Kresoxim-Methyl	0,010	Oxydemeton-Methyl	0,010	Prothiophos	0,010
Flutianil	0,010	Lenacil	0,010	Oxydemeton-Methyl		Pyflubumid	0,010
Flutolanil	0,010	Leptophos	0,010	(Summenparameter)**		Pymetrozin	0,010
Flutriafol	0,010	Linuron	0,010	Oxyfluorfen	0,010	Pyraclostrobin	0,010
Fluvalinat (Summe der Isomere)	0,010	Lufenuron	0,010	Oxymatrin	0,010	Pyraflufen	0,010
Fluxapyroxad	0,010	Malaaxon	0,010	Paclobutrazol (Summe der Isomere)	0,010	Pyraflufen-Ethyl	0,010
FM-6-1	0,010	Malathion	0,010	Paraaxon	0,010	Pyraflufen-Ethyl	
Folpet	0,010	Malathion (Summenparameter)**		Parathion	0,010	(Summenparameter)**	
Folpet (Summenparameter)**		Mandestrobin	0,010	Parathion-Methyl	0,010	Pyrazophos	0,010
Fonofos	0,010	Mandipropamid (Summe der Isomere)	0,010	Parathion-Methyl		Pyrethrin (Summenparameter)**	
Forchlorfenuron	0,010	Matrin	0,010	(Summenparameter)**		Pyrethrin I	0,010
Formetanat	0,010	MCPA (freie Säure)*	0,010	Pebulat	0,010	Pyrethrin II	0,010
Formetanat (als Formetanat-HCl)**		MCPA (Summenparameter)**		Penconazol (Summe der Isomere)	0,010	Pyridaben	0,010
Formothion	0,010	MCPB (freie Säure)*	0,010	Pencycuron	0,010	Pyridafof (CL 9673)	0,010
Fosthiazat	0,010	Mecarbam	0,010	Pencycuron (Summenparameter)**		Pyridalyl	0,010
Fuberidazol	0,010	Mecoprop (Summe der Isomere)	0,010	Pencycuron-PB-amin	0,010	Pyridaphenthion	0,010
Furalaxyl	0,010	Mefenpyr-diethyl	0,010	Pendimethalin	0,010	Pyridat	0,010
Furathiocarb	0,010	Mefentrifluconazol	0,010	Penflufen	0,010	Pyrifeno	0,010
				Pentachloranilin	0,010	Pyrimethanil	0,010
						Pyriofenon	0,010

Pyriproxyfen	0,010	Thionazin	0,010	Zusatztuntersuchungen <i>additional analyses:</i>	Gesamtbromid, BG/RL 1 mg/kg
Quinalphos	0,010	Thiophanat-Methyl	0,010	(nicht in GALAB 500Plus)	Chlorid 5 mg/kg
Quinclorac	0,010	Tiocardbazil	0,010	(not enclosed in GALAB 500Plus)	Bromid-Chlorid-Verhältnis
Quinmerac	0,010	Tolclofos-Methyl	0,010		
Quinmerac (Summenparameter)**		Tolifenpyrad	0,010	Saure Herbizide (alkalische Hydrolyse), <i>phenoxyalkanecarboxylic acids (alkaline hydrolysis),</i>	Glyphosat, AMPA, Glufosinat, N-Acetyl Glufosinat, MPPA LC-MS/MS
Quinoclamim	0,010	Tolylfluanid	0,010		SOP-0657;2022-12
Quinoxifen	0,010	Tolylfluanid (Summenparameter)**		mittels LC-MS/MS	BG/RL je 0,01 mg/kg
Quintozen	0,010	Topramezon	0,005	DIN EN 15662:2018-07 Modulares QuEChERS-	
Quintozen (Summenparameter)**		Tralkoxydime	0,010	Verfahren, alkalische Hydrolyse	
Quizalofof (Summe der Isomere)	0,010	Transfluthrin	0,010	BG/RL je 0,01 mg/kg	Nikotin
(freie Säure)*		Triadimefon	0,010	2,4,5-T	LC-MS/MS
Quizalofof (Summenparameter)**		Triadimenol (Summe der Isomere)	0,010	2,4,5 -T (Summenparameter)**	SOP-196: 2024-06
Quizalofof-Ethyl*	0,010	Triallat	0,010	2,4-D (Summenparameter)**	BG/RL 0,01mg/kg
Rabenzazol	0,010	Triamiphos	0,010	2,4-DB	
Resmethrin (Summe der Isomere)	0,010	Triasulfuron	0,010	2,4-DB (Summenparameter)**	Nitrat
Rimsulfuron	0,010	Triazamat	0,010	2,4-DP (Dichlorprop)	HPLC/IC
Rotenon	0,010	Triazophos	0,010	2,4-DP (Dichlorprop) (Summenparameter)**	Gemäß SOP-0570, Version 2.0; June 2020
Saflufenacil	0,010	Tribenuron-Methyl	0,010	Acibenzolar-S-methyl (Summenparameter)**	BG/RL 1 mg/kg
Saflufenacil (Summenparameter)**		Trichlorfon	0,010	Acibenzolar-Säure	
Saflufenacil, M800H11	0,010	Trichloronat	0,010	Bentazon	Paraquat, Diquat
Saflufenacil, M800H35	0,010	Triclopyr	0,010	Bentazon (Summenparameter)**	LC-MS/MS
Sebutylhazin	0,010	Tricyclazol	0,010	Bentazon-6-OH	SOP-Nr. 495:2022-10; EU-SRM QuPPe 2020-02
Sedaxan	0,010	Tridemorph	0,010	Bentazon-8-OH	BG/RL je 0,01 mg/kg
Sethoxydim	0,010	Trietazin	0,010	Bispyribac	
Siduron	0,010	Trifloxystrobin	0,010	Bispyribac (Summenparameter)**	Perchlorat, Chlorat
Silaneophan	0,010	Triflumizol	0,010	Fluazifop	LC-MS/MS
Silthiofam	0,010	Triflumizol (Summenparameter)**		Fluazifop (Summenparameter)**	SOP-Nr. 495:2022-10; EU-SRM QuPPe 2020-02
Simazin	0,010	Triflururon	0,010	Fluroxypyr (Summenparameter)**	BG/RL je 0,01 mg/kg
Sintofen	0,010	Trifluralin	0,010	Fluroxypyr	
Spinetoram	0,010	Triflursulfuron	0,010	Haloxyfop	Phosphan
Spinosad (Summenparameter)**		Triforin	0,010	Haloxyfop (Summenparameter)**	GC-MSD
Spinosyn A	0,010	Trinexapac	0,010	MCPA	SOP-0559:2019-05
Spinosyn D	0,010	Triticonazol	0,010	MCPA (Summenparameter)**	BG/RL 0,01 mg/kg
Spirodiclofen	0,010	Tritosulfuron	0,010	MCPB	
Spiromesifen	0,010	Uniconazol	0,010	Phenylphenol, ortho-	Pyridat, Pyridafol, Pyridat (Summenparameter)**
Spirotetramat	0,010	Valifenalat	0,010	Phenylphenol, ortho- (Summenparameter)**	nach Hydrolyse
Spirotetramat		Vamidothion	0,010	Quizalofof (Summe der Isomeren)	LC-MS/MS
(Summenparameter)**		Vinclozolin	0,010	Quizalofof (Summenparameter)**	DIN EN 15662:2018-07 Modulares QuEChERS-
Spirotetramat-enol	0,010	Zoxamid	0,010		Verfahren, alkalische Hydrolyse
Spirotetramat-enol-Glc	0,010				BG/RL 0,01mg/kg
Spirotetramat-ketohydroxy	0,010			Chlormequat, Mepiquat	
Spirotetramat-monohydroxy	0,010			LC-MS/MS	
Spiroxamin (Summe der Isomere)	0,010			SOP-Nr. 495:2022-10; EU-SRM QuPPe 2020-02	Trimethylsulfonium (Trimesium)
Sulcotrion	0,010			BG/RL je 0,01 mg/kg	LC-MS/MS
Sulfentrazon	0,010			Chlormequat, BG/RL 0,01 mg/kg	SOP-Nr. 495:2022-10; EU-SRM QuPPe 2020-02
Sulfosulfuron	0,010			Mepiquat, BG/RL 0,01 mg/kg	BG/RL 0,01mg/kg
Sulfotep	0,010			Chlormequat (Summenparameter)**,	
Sulfoxaflor	0,010			Mepiquat (Summenparameter)**	
Sulprofos	0,010				
Swep	0,010			Difluoressigsäure, Trifluoressigsäure (DFA, TFA)	
Tebuconazol	0,010			LC-MS/MS, SOP-Nr. 495:2022-10; EU-SRM QuPPe	
Tebufenozid	0,010			2020-02	
Tebufenpyrad	0,010			BG/RL je 0,01 mg/kg	
Tebupirimphos	0,010				
Tebutam	0,010			Dithiocarbamate, berechnet als CS ₂ /calculated as CS ₂	
Tecnazen	0,010			gem. ASU §64 LFGB L 00.00-49/2 modif.: GC-MSD	
Teflubenzuron	0,010			BG/RL 0,01 mg/kg,	
Tefluthrin (Summe der Isomere)	0,010				
Tembotrion	0,010			Ethephon	
Temephos	0,010			LC-MS/MS	
Tepraloxymid	0,010			SOP-Nr. 495:2022-10; EU-SRM QuPPe 2020-02	
Terbacil	0,010			BG/RL 0,01 mg/kg	
Terbufos	0,010				Ethylenoxid (Summenparameter)**
Terbumeton	0,010			GC-MSMS, SOP-0653:02-2022	
Terbutylhazin	0,010			BG/RL 0,01mg/kg	
Terbutryn	0,010				Ethylenoxid, 2-Chlorethanol, Ethylenoxid
Tetrachlorvinphos	0,010				(Summe)
Tetraconazol (Summe der Isomere)	0,010			GC-MSMS, SOP-0653:02-2022	
Tetradifon	0,010			BG/RL 0,01mg/kg	
Tetrahydrophthalimid	0,010			Ethylenoxid (Summenparameter)**	
Tetramethrin	0,010			Ethylenoxid, BG/RL 0,01 mg/kg	
Tetrasul	0,010			2 Chlorethanol, BG/RL 0,01 mg/kg	
TFNA	0,010				
TFNG	0,010			Fosetyl (Summenparameter)**	
Thiabendazol	0,010			LC-MS/MS	
Thiacloprid	0,010			SOP-Nr. 495:2022-10; EU-SRM QuPPe 2020-02	
Thiamethoxam	0,010			Fosetyl, BG/RL 0,01 mg/kg	
Thidiazuron	0,010			Phosphonsäure, BG/RL 0,01 mg/kg	
Thifensulfuron-Methyl	0,010				
Thiodicarb	0,010				Gesamtbromid
Thiofanox	0,010				GC-MSD
Thiofanox-sulfon	0,010				ASU L 00.00-36/1:2004-07,
Thiofanox-sulfoxid	0,010				BG/RL 1 mg/kg
Thiometon	0,010				
					Gesamtbromid, Chlorid, Br/Cl-Verhältnis
					GC-MSD
					ASU L 00.00-36/1:2004-07,