

Substanzname	Berichtsgrenze/ reporting limit (RL) [mg/kg]						
1,2,3,4-Tetrachlorbenzol	0,010	Bifenox	0,010	Clofentezin	0,050	Dinoseb	0,10
1,2,4,5-Tetrachlorbenzol	0,010	Bifenthrin (Summe der Isomere)	0,010	Clomazon	0,010	Dinoseb (Summenparameter)**	
1,4-Dimethylnaphthalin	0,010	Binapacryl	0,010	Clopyralid	0,010	Dinoseb-Acetat	0,010
1-Naphthylacetamid (1-NAD)	0,010	Biphenyl	0,010	Cloquintocet-mexyl	0,010	Dinotefuran	0,010
1-Naphthyllessigsäure	0,010	Bispyribac (freie Säure)*	0,010	Clothianidin	0,010	Dioxacarb	0,010
1-Naphthyllessigsäure	0,010	Bispyribac		Coumaphos	0,010	Dioxathion	0,010
(Summenparameter)**		(Summenparameter)**		Crimidin	0,010	Diphenamid	0,010
2,3,5-Trimethacarb/3,4,5-Trimethacarb	0,010	Bittertanol (Summe der Isomere)	0,010	Crufomat	0,010	Diphenylamin	0,010
2,4-D (freie Säure)*	0,010	Boscalid	0,010	Cyanazin	0,010	Dipropetryn	0,010
2,4-D (Summenparameter)**		Bromfenvinphos	0,010	Cyanophos	0,010	Disulfoton	0,010
2,4-DB (freie Säure)*	0,010	Bromocyclen	0,010	Cyantraniliprol	0,010	Disulfoton	
2,4-DB (Summenparameter)**		Bromophos-Ethyl	0,010	Cyazofamid	0,010	(Summenparameter)**	
2,4-DP (freie Säure)*	0,010	Bromophos-Methyl	0,010	Cyclanilid	0,010	Disulfoton-Sulfon	0,010
2,4-DP (Summenparameter)**		Bromoxynil	0,050	Cycloat	0,010	Disulfoton-Sulfoxid	0,010
3-Decen-2-On	0,010	Bromopropylat	0,010	Cyflufenamid (Summe der Isomere)	0,010	Dithianon	0,010
479M04	0,010	Bromuconazol	0,010	Cyflumetofen (Summe der Isomere)	0,010	Diuron	0,010
479M08	0,010	(Diastereoisomere)		Cyfluthrin (Summe aller Isomere)	0,010	DNOC	0,010
479M16	0,010	BTS 44595	0,010	Cyhalothrin, lambda-/gamma-		Dodemorph	0,010
4-Bromphenylharnstoff	0,010	BTS 44596	0,010	Cyhexatin (Summe Azocyclotin und Cyhexatin)	0,010	Dodin	0,050
4-CPA	0,010	Bupirimat	0,010	Cymiazol	0,010	Edifenphos	0,010
8-Hydroxychinolin	0,010	Buprofezin	0,010	Cymoxanil	0,050	Emamectin B1a (freie Base)	0,002
Abamectin (Summenparameter)**		Butafenacil	0,010	Cypermethrin (Summe aller Isomere)	0,010	Endosulfan	
Acephat	0,050	Cadusafos	0,010	DDAC (Summenparameter)**		Endosulfan, alpha-	0,010
Acequinocyl	0,010	Captan	0,010	DDAC 10	0,010	Endosulfan, beta-	0,010
Acetamiprid	0,010	Captan (Summenparameter)**		DDAC 12	0,010	Endosulfansulfat	0,010
Acetochlor	0,010	Carbaryl	0,010	DDAC 8	0,010	Epoxiconazol	0,010
Acibenzolar-Säure (freie Säure)*	0,010	Carbendazim	0,010	DDD, p,p'-	0,010	EPTC	0,050
Acibenzolar-S-methyl	0,010	Carbendazim (Summenparameter)**		DDD, p,p'-	0,010	Ethiofencarb	0,010
Acibenzolar-S-methyl (Summenparameter)**		Carboxin (Summenparameter)**		DDT (Summenparameter)**		Ethion	0,010
Acionifen	0,010	Carboxin-Sulfoxid	0,010	DDT, o,p'-	0,010	Ethiprol	0,010
Acrinathrin	0,010	Carfentrazon-Ethyl	0,010	DEF	0,010	Ethofumesat	0,010
Alachlor	0,010	Cetrimoniumchlorid	0,010	Deltamethrin	0,010	Ethofumesat-2-keto	0,010
Aldicarb	0,010	Chinomethionat	0,010	Denatoniumbenzoat	0,010	Ethofumesat-Carbonsäure	0,010
Aldicarb (Summenparameter)**		Chlorantraniliprol	0,010	Desethyl-Simazin	0,10	Ethoprophos	0,010
Aldicarb-sulfoxid	0,010	Chlorbensid	0,10	Desmedipham	0,10	Ethoxyquin	0,010
Aldoxycarb	0,010	Chlorbenzilat	0,010	Desmetryn	0,010	Etofenprox	0,010
Aldrin	0,010	Chlorbufam	0,010	Diaphenthiuron	0,010	Etozaxol	0,050
Allethrin	0,010	Chlordan (cis-/trans-)**		Diallat	0,050	Etrifos	0,010
Ametryn	0,010	Chlordan, cis-	0,010	Diazinon	0,010	Famoxadon	0,050
Amidosulfuron	0,010	Chlordan, trans-	0,010	Dicamba	0,010	Famphur	0,010
Aminocarb	0,010	Chlordecon	0,010	Dichlobenil	0,010	Fenamidon	0,010
Amisulbrom	0,010	Chlorfenapyr	0,010	Dichlofenthion	0,010	Fenamiphos	0,010
Ancymidol	0,010	Chlorfenprop-Methyl	0,010	Dichlofluanid	0,010	(Summenparameter)**	
Anthrachinon	0,010	Chlorfenson	0,010	Diclobutrazol	0,010	Fenazaquin	0,010
Atrazin	0,010	Chlorfenvinphos	0,010	Diclofop	0,010	Fenbuconazol	0,010
Avermectin B1a	0,006	Chloridazon	0,010	Diclofop-Methyl	0,010	Fenbutatinoxid	0,050
Avermectin B1b	0,006	(Summenparameter)**		Diclofop-Methyl (Summenparameter)**		Fenchlorphos	0,010
Azadirachtin	0,010	Chlormephos	0,010	Dicloran	0,010	Fenchlorphos-oxon	0,10
Azamethiphos	0,010	Chloroneb	0,010	Dicofol (Summe aus p,p- und o,p-Isomere)	0,010	Fenfluthrin	0,010
Azinphos-Ethyl	0,010	Chloroxuron	0,010	Dieldrin	0,010	Fenhexamid	0,010
Azinphos-Methyl	0,010	Chlorpropham	0,010	Diethofencarb	0,010	Fenitrothion	0,010
Aziprotryn	0,010	Chlorpyrifos	0,010	Diethyltoluamid (DEET)	0,010	Fenobucarb	0,010
Azoxystrobin	0,010	Chlorpyrifos-Methyl	0,010	Difenoconazol	0,010	Fenoprop (2,4,5-TP)	0,010
BAC (Summenparameter)**		Chlorsulfuron	0,010	Difenoxuron	0,010	Fenothiocarb	0,010
BAC 10	0,010	Chlorthal-Dimethyl	0,010	Diffubenzuron	0,010	Fenoxaprop-Ethyl	0,10
BAC 12	0,010	Chlorthalonil	0,010	Diffufenican	0,010	Fenoxycarb	0,010
BAC 14	0,010	Chlorthiamid	0,010	Dimefox	0,010	Fenpicoxamid	0,010
BAC 16	0,010	Chlorthion	0,010	Dimethachlor	0,010	Fenpropathrin	0,010
BAC 18	0,010	Chlortoluron	0,050	Dimethenamid (Summe aller Isomere)	0,010	Fenpropidin	0,010
BAC 8	0,010	Chlozolinat	0,050	Dimethipin	0,010	Fenpropimorph (Summe der Isomere)	0,010
Benalaxyl (Summe der Isomere)	0,010	Cinerin I	0,010	Dimethoat	0,010	Fenpyroximat	0,010
Bendiocarb	0,010	Cinerin II	0,010	Dimethomorph (Summe der Isomere)	0,010	Fenson	0,010
Benfluralin	0,010	Cinosulfuron	0,010	Dimethylaminosulfotoluidid	0,010	Fensulfothion	0,010
Benfuracarb	0,010	Clethodim	0,010	Dimetilan	0,010	Fensulfothion-sulfon	0,010
Bensulfuron-Methyl	0,010	(Summenparameter)**		Dimoxystrobin	0,010	Fenthion (Summenparameter)**	
Bentazon	0,050	Clethodim sulfon	0,010	Diniconazol (Summe der Isomere)	0,010	Fenthion-Oxon	0,010
Bentazon (Summenparameter)**		Clethodim sulfoxid	0,010	Dinitramin	0,010	Fenthion-Oxonsulfon	0,010
Bentazon-6-OH*	0,010	Clodinafop				Fenthion-Oxonsulfoxid	0,010
Bentazon-8-OH*	0,010	(Summenparameter)**				Fenthion-sulfon	0,010
Benthiavalicarb-Isopropyl	0,010	Clodinafop-Propargyl	0,010			Fentin	0,010
Benzovindiflupyr	0,010					Fenuron	0,010
Benzoylprop-Ethyl	0,010					Fenvalerat/Esfenvalerat (Summe)	0,010
BH 518-2	0,010						
BH 518-4	0,010						
Bifenazat	0,010						
Bifenazat (Summenparameter)**							
Bifenazat-diazen	0,010						

Fipronil	0,005	Imazethapyr	0,010	Mevinphos (Summe der E- und Z-Isomere)	0,010	Prochloraz	0,010
Fipronil (Summenparameter)**		Imidacloprid	0,010	Milbectin		Prochloraz (Summenparameter)**	
Fipronil-Sulfoxid	0,005	Indaziflam	0,010	(Summenparameter)**		Procymidon	0,010
Fipronil-sulfon	0,005	Indoxacarb (Summe der Isomere)	0,010	Milbectin A3	0,010	Profenofos	0,010
Flazasulfuron	0,010	Iodofenphos	0,010	Milbectin A4	0,010	Profluralin	0,010
Flometoquin	0,010	Iodosulfuron-Methyl	0,010	Mirex	0,010	Profoxydim	0,010
Flonicamid	0,010	Ioxynil	0,050	Molinat	0,010	Promecarb	0,010
Flonicamid (Summenparameter)**		Iprobenfos	0,010	Monocrotophos	0,050	Prometon	0,010
Florpyrauxifen-benzyl	0,010	Iprodion	0,010	Monolinuron	0,010	Prometryn	0,010
Fluacrypyrim	0,010	Iprovalicarb	0,010	Monuron	0,010	Propachlor	0,010
Fluazifop (freie Säure)*	0,010	Isobenzan	0,010	Myclobutanil	0,010	(Summenparameter)**	
Fluazifop (Summenparameter)**		Isocarboxin	0,010	Napropamid	0,010	Propachlor	
Fluazifop-Butyl*	0,010	Isodrin	0,010	Neburon	0,010	Propamocarb	0,050
Fluazinam	0,050	Isofetamid	0,010	Nicosulfuron	0,010	Propaquizafop*	0,010
Flubendiamid	0,010	Isomethiozin	0,010	Nitralin	0,010	Propargit	0,010
Fluchloralin	0,010	Isoprocab	0,010	Nitrapyrin	0,010	Propazin	0,010
Flucythrinate	0,010	Isopropalin	0,010	Norflurazon	0,010	Propetamphos	0,010
Fludioxonil	0,010	Isoprothiolan	0,010	Novaluron (Summe der Isomere)	0,010	Propham	0,010
Fluensulfon	0,010	Isoproturon	0,010	Ofurac	0,010	Propiconazol	0,010
Flufenacet	0,010	Isoprazam	0,010	Omethoat	0,010	Propoxur	0,005
Flufenacet (Summenparameter)**		Isoxaben	0,010	Oxadiazon	0,050	Propyzamid	0,010
Flufenacet oxalat	0,010	Isoxadifen-Ethyl	0,10	Oxadixyl	0,010	Proquinazid	0,010
Flufenacet Sulfonsäure	0,010	Isoxaflutol	0,010	Oxamyl	0,001	Prosulfocarb	0,010
Flufenacet thioglycolate sulfoxid	0,010	Isoxaflutol (Summenparameter)**		Oxathiapropilin	0,010	Prosulfuron	0,010
Flufenoxuron	0,010	Isoxaflutol diketonitril	0,010	Oxycarboxin	0,010	Prothioconazol	0,010
Flumethrin	0,010	Isoxathion	0,010	Oxydemeton-Methyl	0,050	Prothioconazol-Desthio	0,010
Flumioxazin	0,10	Jasmodin I	0,010	(Summenparameter)**		(Summenparameter)**	
Flumeturon	0,005	Jasmodin II	0,010	Oxyfluorfen	0,050	Prothiophos	0,010
Flupicolid	0,010	Kresoxim-Methyl	0,010	Oxymatrin	0,010	Pyflubumid	0,010
Fluopyram	0,010	Lenacil	0,10	Paclobutrazol (Summe der Isomere)	0,010	Pymetrozin	0,010
Fluoxastrobin	0,010	Leptophos	0,010	Paraaxon	0,010	Pyraclostrobin	0,010
Flupyradifuron	0,010	Linuron	0,010	Paraaxon-Methyl	0,010	Pyraflufen	0,010
Fluquinconazol	0,010	Lufenuron	0,010	Parathion-Methyl	0,010	Pyraflufen-Ethyl	0,010
Flurochloridone	0,010	Malaoxon	0,010	Parathion-Methyl (Summenparameter)**		(Summenparameter)**	
Fluroxypyr (freie Säure)*	0,010	Malathion	0,010	Pebulat	0,010	Pyrazophos	0,010
Fluroxypyr (Summenparameter)**		Malathion (Summenparameter)**		Penconazol (Summe der Isomere)	0,010	Pyrethrin (Summenparameter)**	
Flusilazol	0,010	Mandestrobin	0,010	Pencycuron	0,010	Pyrethrin I	0,010
Flutianil	0,010	Mandipropamid (Summe der Isomere)	0,010	Pencycuron (Summenparameter)**		Pyrethrin II	0,010
Flutriafol	0,010	Matrin	0,010	Pencycuron-PB-amin	0,010	Pyradaben	0,010
Fluvalinat (Summe der Isomere)	0,010	MCPA (freie Säure)*	0,010	Pendimethalin	0,010	Pyradafol (CL 9673)	0,010
Fluxapyroxad	0,010	MCPB (freie Säure)*	0,010	Penflufen	0,010	Pyradaphenthion	0,010
FM-6-1	0,010	Mecarbam	0,010	Pentachloranilin	0,010	Pyradat	0,010
Folpet	0,010	Mecoprop (Summe der Isomere)	0,10	Pentachloranisole	0,010	Pyrifenox	0,010
Folpet (Summenparameter)**		Mefenpyr-diethyl	0,010	Pentachlorbenzol	0,010	Pyrimethanil	0,010
Fonofos	0,010	Mefentrifluconazol	0,010	Pentachlorphenol	0,010	Pyrifonon	0,010
Forchlorfenuron	0,010	Mepanipyrim	0,010	Pentachlorphenol	0,010	Pyriproxyfen	0,010
Fosthiazat	0,010	Mephosfolan	0,010	Permethrin (Summe der Isomere)	0,010	Quinalphos	0,010
Fuberidazol	0,010	Merphos	0,010	Perthan	0,010	Quinclorac	0,010
Furalaxyl	0,010	Metalaxyl (Summe der Isomere)	0,010	Pethoxamid	0,010	Quinmerac	0,10
Furathiocarb	0,010	Metaldehyd	0,010	Phenmedipham	0,010	(Summenparameter)**	
Furmecycloz	0,010	Metazachlor	0,010	Phenothrin	0,050	Quinoclamid	0,010
Genite	0,010	(Summenparameter)**		Phenthoat	0,010	Quinoxifen	0,010
Halfenprox	0,010	Metconazol (Summe der Isomere)	0,010	Phenylphenol, ortho-	0,010	Quintozen	0,010
Halosulfuron-Methyl	0,010	Methabenzthiazuron	0,050	Phenylphenol, ortho- (Summenparameter)**		Quintozen	
Haloxypol (freie Säure)*	0,010	Methacryfos	0,010	Phorac (Summenparameter)**		(Summenparameter)**	
Haloxypol (Summenparameter)**		Methidathion	0,010	Phorac-Oxon	0,010	Quizalofop (Summe der Isomere) (freie Säure)*	0,010
Haloxypol-ethoxyethylster*	0,010	Methiocarb	0,010	Phorac-Oxon-Sulfon	0,010	Quizalofop	
Haloxypol-methylester*	0,010	Methiocarb-sulfon	0,010	Phosalan	0,010	(Summenparameter)**	
HCH, alpha-	0,010	Methiocarb-sulfoxid	0,010	Phosfolan	0,010	Quizalofop-Ethyl*	0,010
HCH, beta-	0,010	Methomyl	0,010	Phosmet	0,005	Rabenzazol	0,010
HCH, delta-	0,010	Methoprotin	0,010	Phosphamidon	0,010	Resmethrin (Summe der Isomere)	0,010
HCH, epsilon-	0,010	Methoxychlor	0,10	Phoxim	0,010	Rimsulfuron	0,050
HCH, gamma-	0,010	Methoxyfenozid	0,010	Phthalimid	0,010	Rotenon	0,010
Heptachlor	0,010	Metobromuron	0,010	Picloram	0,010	Saflufenacil	0,010
Heptachlor (Summenparameter)**		(Summenparameter)**		Picolinafen	0,010	Saflufenacil	
Heptachlorepoxyd, cis-	0,010	Metolachlor (Summe der Isomere)	0,010	Picoxystrobin	0,010	(Summenparameter)**	
Heptachlorepoxyd, trans-	0,010	Mitolcarb	0,010	Piperonylbutoxid	0,010	Saflufenacil, M800H11	0,010
Heptenophos	0,010	Metominostrobin	0,010	Piperophos	0,010	Saflufenacil, M800H35	0,010
Hexachlorbenzol	0,010	Metoxuron	0,010	Pirimicarb	0,010	Sebutylazin	0,010
Hexaconazol	0,010	Metsulfuron-Methyl	0,050	Pirimiphos-Ethyl	0,010	Sedaxan	0,010
Hexazinon	0,010			Pirimiphos-Methyl	0,010	Sethoxydim	0,010
Hexythiazox	0,010			Prallethrin	0,010	Siduron	0,010
Hymexazol	0,010			Pretilachlor	0,010	Silthiofamid	0,010
Icaridin	0,010					Simazin	0,050
Imazalil	0,050					Sintofen	0,010
Imazamox	0,010					Spinetoram	0,010
Imazapyr	0,010					Spinosad (Summenparameter)**	
Imazaquin	0,010					Spinosyn A	0,010

Spinosyn D	0,010	*Für die Bewertung des	Chlormequat, Mepiquat
Spirodiclofen	0,010	Rückstandshöchstgehaltes ist die	LC-MS/MS
Spiromesifen	0,010	Zusatzuntersuchung der sauren Herbizide	SOP-Nr. 495:2022-10; EU-SRM QuPpE 2020-02
Spirotetramat	0,010	(alkalische Hydrolyse) notwendig! <i>For the legal</i>	BG/RL je 0,01 mg/kg
Spirotetramat		<i>assessment of the MRL an additional test of</i>	Chlormequat, BG/RL 0,01 mg/kg
(Summenparameter)**		<i>phenoxyalkancarboxylic acids (alkaline</i>	Mepiquat, BG/RL 0,01 mg/kg
Spirotetramat-enol	0,010	<i>hydrolysis) is necessary.</i>	Chlormequat (Summenparameter)**,
Spirotetramat-enol-Glc	0,010		Mepiquat (Summenparameter)**
Spirotetramat-ketohydroxy	0,010	**berechnet/calculated	
Spirotetramat-monohydroxy	0,010	Methodik Method	Difluoressigsäure, Trifluoressigsäure (DFA, TFA)
Spiroxamin (Summe der	0,010	DIN EN 15662:2018-07 Modulares QuEChERS-	LC-MS/MS, SOP-Nr. 495:2022-10; EU-SRM QuPpE
Isomere)		Verfahren	2020-02
Sulcotrion	0,010		BG/RL je 0,01 mg/kg
Sulfentazon	0,010	Messsystem Technical equipment	Dithiocarbamate, berechnet als CS ₂ /calculated as
Sulfosulfuron	0,010	LC-MS/MS	CS ₂
Sulfotep	0,010	GC-MSD/GC-MS/MS	gem. ASU §64 LFGB L 00.00-49/2 modif.: GC-MSD
Sulfoxaflor	0,010	GC-NCI	BG/RL 0,01 mg/kg,
Sulprofos	0,010		
Tebuconazol	0,010		Ethephon
Tebufenozid	0,010		LC-MS/MS
Tebufenpyrad	0,010	Zusatzuntersuchungen additional analyses:	SOP-Nr. 495:2022-10; EU-SRM QuPpE 2020-02
Tebupirimphos	0,050	(nicht in GALAB 500Plus komplexe Matrices	BG/RL 0,01 mg/kg
Tebutam	0,010	enthalten)	
Teflubenzuron	0,010	(not enclosed in GALAB 500Plus complex	Ethylenoxid (Summenparameter)**
Tefluthrin (Summe der Isomere)	0,010	Matrices)	GC-MSMS, SOP-0653:02-2022
Tepraloxymid	0,010	Saure Herbizide (alkalische Hydrolyse),	BG/RL 0,01mg/kg
Terbufos	0,010	phenoxyalkancarboxylic acids (alkaline	Ethylenoxid, 2-Chlorethanol, Ethylenoxid
Terbumeton	0,010	hydrolysis),	(Summe)
Terbuthylazin	0,010	mittels LC-MS/MS	GC-MSMS, SOP-0653:02-2022
Terbutryn	0,010	BG/RL je 0,01 mg/kg	BG/RL 0,01mg/kg
Tetrachlorvinphos	0,010	2,4,5-T	Ethylenoxid (Summenparameter)**
Tetraconazol (Summe der	0,010	2,4,5-T (Summenparameter)**	Ethylenoxid, BG/RL 0,01 mg/kg
Isomere)		2,4-D (Summenparameter)**	2 Chlorethanol, BG/RL 0,01 mg/kg
Tetradifon	0,010	2,4-D	
Tetrahydrophthalimid	0,010	2,4-DB	Fosetyl (Summenparameter)**
Tetramethrin	0,010	2,4-DB (Summenparameter)**	LC-MS/MS
TFNA	0,010	2,4-DP (Dichlorprop)	SOP-Nr. 495:2022-10; EU-SRM QuPpE 2020-02
TFNG	0,010	2,4-DP (Dichlorprop) (Summenparameter)**	Fosetyl, BG/RL 0,01 mg/kg
Thiabendazol	0,050	Acibenzolar-S-methyl (Summenparameter)**	Phosphonsäure, BG/RL 0,01 mg/kg
Thiacloprid	0,010	Acibenzolar-Säure	
Thiamethoxam	0,050	Bentazon	Gesamtbromid
Thidiazuron	0,010	Bentazon (Summenparameter)**	GC-MSD
Thifensulfuron-Methyl	0,050	Bentazon-6-OH	ASU L 00.00-36/1:2004-07,
Thiodicarb	0,010	Bentazon-8-OH	BG/RL 1 mg/kg
Thiophanat-Methyl	0,10	Bispyribac	
Tiocarbazil	0,010	Bispyribac (Summenparameter)**	Gesamtbromid, Chlorid, Br/Cl-Verhältnis
Tolclofos-Methyl	0,010	Fluazifop	GC-MSD
Tolfenpyrad	0,010	Fluazifop (Summenparameter)**	ASU L 00.00-36/1:2004-07,
Tolyfluanid	0,010	Fluroxyppyr	Gesamtbromid, BG/RL 1 mg/kg
(Summenparameter)**		Fluroxyppyr (Summenparameter)**	Chlorid 5 mg/kg
Topramezon	0,005	Haloxypfop	Bromid-Chlorid-Verhältnis
Tralkoxydime	0,010	Haloxypfop (Summenparameter)**	
Transfluthrin	0,010	MCPA	Glyphosat
Triadimefon	0,010	MCPA (Summenparameter)**	LC-MS/MS
Triadimenol (Summe der	0,010	MCPB	SOP-0232:2011-06
Isomere)		Phenylphenol, ortho-	BG/RL je 0,01 mg/kg
Triallat	0,010	Phenylphenol, ortho- (Summenparameter)**	
Triamiphos	0,010	Quiazalofop (Summe der Isomeren)	Nitrat
Triasulfuron	0,050	Quiazalofop (Summenparameter)**	HPLC/IC
Triazamat	0,010		Gemäß SOP-0570, Version 2.0; June 2020
Triazophos	0,010		BG/RL 1 mg/kg
Tribenuron-Methyl	0,010		
Trichlorfon	0,050		Paraquat, Diquat
Trichloronat	0,010		mittels LC-MS/MS
Triclopyr	0,010		SOP-Nr. 495:2022-10; EU-SRM QuPpE 2020-02
Tricyclazol	0,010		BG/RL je 0,01 mg/kg
Tridemorph	0,010		
Trietazin	0,010		Perchlorat, Chlorat
Trifloxystrobin	0,010		LC-MS/MS
Triflumizol	0,010		SOP-Nr. 495:2022-10; EU-SRM QuPpE 2020-02
Triflumizol			BG/RL je 0,01 mg/kg
(Summenparameter)**			
Triflumuron	0,050		Phosphan
Trifluralin	0,010		GC-MSD
Triflusulfuron	0,010		SOP-0559:2019-05
Triforin	0,010		BG/RL 0,01 mg/kg
Trinexapac	0,010		
Tritosulfuron	0,010		Pyridat, Pyridafol, Pyridat (Summenparameter)**
Vamidothion	0,010		nach Hydrolyse
Vinclozolin	0,010		LC-MS/MS
Zoxamid	0,010		DIN EN 15662:2018-07 Modulares QuEChERS-
			Verfahren, alkalische Hydrolyse
			BG/RL 0,01mg/kg
			Sulfit
			Gem. § 64 LFGB ASU L 00.00 46/1, November 1999
			BG/RL 10 mg/kg
			Trimethylsulfonium (Trimesium)
			LC-MS/MS
			SOP-Nr. 495:2022-10; EU-SRM QuPpE 2020-02
			BG/RL 0,01mg/kg